

4.1.08-00 Peltier heat pump



What you can learn about ...

- Peltier effect
- Heat pipe
- Thermoelectric e. m. f.
- Peltier coefficient
- Cooling capacity
- Heating capacity
- Efficiency rating
- Thomson coefficient
- Seebeck coefficient
- Thomson equations
- Heat conduction
- Convection
- Forced cooling
- Joule effect


Principle:

The cooling capacity heating capacity and efficiency rating of a Peltier heat pump are determined under different operating conditions.

What you need:

Thermogenerator	04366.00	1
Flow-through heat exchanger	04366.01	1
Air cooler	04366.02	1
Heating coil with sockets	04450.00	1
Distributor	06024.00	1
Rheostats, 33 Ω , 3.1 A	06112.02	1
Connecting plug, pack of 2	07278.05	1
Power supply, universal	13500.93	1
Digital multimeter 2010	07128.00	4
Stopwatch, digital, 1/100 s	03071.01	1
Hot/cold air blower, 1700 W	04030.93	1
Laboratory thermometers, -10...+100°C	38056.00	1
Precision mercury thermometers, -10...+ 50°C	38033.00	2
Rubber tubing, $d_i = 6$ mm, $l = 1$ m	39282.00	1
Universal clamp	37718.00	1
Tripod base -PASS-	02002.55	1
Support rod -PASS-, square, $l = 250$ mm	02025.55	1
Right angle clamp -PASS-	02040.55	1
Connecting cable, 4 mm plug, 32 A, red, $l = 25$ cm	07360.01	3
Connecting cable, 4 mm plug, 32 A, red, $l = 50$ cm	07361.01	3
Connecting cable, 4 mm plug, 32 A, blue, $l = 50$ cm	07361.04	2
Connecting cable, 4 mm plug, 32 A, blue, $l = 75$ cm	07362.04	2
Connecting cable, 4 mm plug, 32 A, red, $l = 75$ cm	07362.01	1
Heat conductive paste, 50 g	03747.00	1

Complete Equipment Set, Manual on CD-ROM included
Peltier heat pump P2410800



Pump cooling capacity as a function of the operating current.

Tasks:

- To determine the cooling capacity P_c the pump as a function of the current and to calculate the efficiency rating η_c at maximum output.
- To determine the heating capacity P_w of the pump and its efficiency rating η_w at constant current and constant temperature on the cold side.
- To determine P_w , η_w and P_c , η_c from the relationship between temperature and time on the hot and cold sides.
- To investigate the temperature behaviour when the pump is used for cooling, with the hot side air-cooled.