

Magnetostriction with the Michelson interferometer 4.3.08-00

NEW



What you can learn about ...

- Interference
- Wavelength
- Diffraction index
- Speed of light
- Phase
- Virtual light source
- Ferromagnetic material
- Weiss molecular magnetic fields
- Spin-orbit coupling

Principle:

With the aid of two mirrors in a Michelson arrangement, light is brought to interference. Due to the magnetostrictive effect, one of the mirrors is shifted by variation in the magnetic field applied to a sample, and the change in the interference pattern is observed.


What you need:

Optical base plate with rubberfeet	08700.00	1
Laser, He-Ne 0.2/1.0 mW, 220 V AC*	08180.93	1
Adjusting support 35 x 35 mm	08711.00	3
Surface mirror 30 x 30 mm	08711.01	4
Magnetic foot for optical base plate	08710.00	7
Holder for diaphragm/ beam splitter	08719.00	1
Beam splitter 1/1, non polarizing	08741.00	1
Lens, mounted, $f = +20$ mm	08018.01	1
Lens holder for optical base plate	08723.00	1
Screen, white, 150 x 150 mm	09826.00	1
Faraday modulator for optical base plate	08733.00	1
Rods for magnetotrixtion, set of 3	08733.01	1
Power supply, universal	13500.93	1
Digital multimeter 2010	07128.00	1
Flat cell battery, 9 V	07496.10	1
Connecting cable, 4 mm plug, 32 A, blue, $l = 50$ cm	07361.04	1

*Alternative:

He/Ne Laser, 5mW with holder	08701.00	1
Power supply for laser head 5 mW	08702.93	1

Complete Equipment Set, Manual on CD-ROM included
Magnetostriction with the Michelson
interferometer **P2430800**



Measuring results of the magnetostriction of nickel with the relative change in length $\Delta l/l$ plotted against applied field strength H .

Tasks:

1. Construction of a Michelson interferometer using separate optical components.
2. Testing various ferromagnetic materials (iron and nickel) as well as a non-ferromagnetic material (copper), with regard to their magnetostrictive properties.