

## 1.5.19-15 Interference of two identical ultrasonic transmitters



## What you can learn about ...

- Longitudinal waves
- Sound pressure
- Huygens' principle
- Interference
- Fraunhofer and Fresnel diffraction


## Principle:

Ultrasonic waves of the same frequency, amplitude and direction of propagation are generated by two sources positioned parallel to each other. The sources can vibrate both in-phase and out-of phase. The angular distribution of the intensity of the waves, which interfere with each other, is automatically recorded using a motor-driven ultrasonic detector and a PC.

## What you need:

Goniometer with reflecting mirror	13903.00	1
Goniometer Operation Unit	13903.99	1
Ultrasound operation unit	13900.00	1
Power supply 5 VDC/2.4 A with DC-socket 2.1 mm	13900.99	1
Ultrasonic transmitter	13901.00	2
Ultrasonic receiver on stem	13902.00	1
Barrel base -PASS-	02006.55	2
Data cable 2 x SUB-D, plug/socket, 9 pole	14602.00	1
Measuring tape, $l = 2$ m	09936.00	1
Screened cable, BNC, $l = 750$ mm	07542.11	1
Adapter BNC socket/4 mm plug pair	07542.27	1
Software Goniometer	14523.61	1
PC, Windows® 95 or higher		

**Complete Equipment Set, Manual on CD-ROM included**  
**Interference of two identical ultrasonic transmitters**  
**P2151915**



Angular distribution of the intensity of two interfering ultrasonic waves having the same phase, amplitude, frequency and direction of propagation.

## Tasks:

1. Determine the angular distribution of the sound pressure of two ultrasonic transmitters vibrating in-phase.
2. Determine the angular positions of the interference minima and compare the values found with those theoretically expected.
3. Repeat the measurements with the two ultrasonic transmitters vibrating out-of-phase.
4. Repeat the first measurement and additionally determine with the angular distribution of the sound pressure of each single transmitter.